




VISIO-TEK™

Wireless digital microscope

THE PORTABLE
DIGITAL MICROSCOPE
WITH INTEGRATED WI-FI

VISIO-MDP5: Dermoscope, Polarized, 1.3 Mpx, 5x-200x, 8 white LEDs
Wi-Fi mode: 1.3 Mpx, USB Mode: 5 Mpx

VISIO-MDP5 is a digital Dermoscope with a polarizer lens which gives you a darker, more saturated, sharper image with increased contrast. Polarizing the light increases and maximizes vision quality.



No Polarized Light



Polarized Light

(FDA LISTED)

FEATURES

- Photo & video Wi-Fi mode resolution: 1.3MPx -1280x1024, 640x480, 320x240
- Photo & video USB mode resolution: 2592x1944 px, 2048x1536 px, 1600x1200 px, 1280x960 px, 1280x1024 px, 640x480 px, 320x240 px.
- Adjustable magnification: 10x ~ 200x
- Frame rate: up to 30fps
- Polarizing Filter
- IEEE 802.11 b/g/n standards compliant
- Wireless transmission distance: up to 5m
- Frequency range: 2.400GHz ~ 2.4835GHz
- Li-ion battery: continuous working time - appx. 2 hours
full charging time - appx. 2 hours
- Camera sensor: CMOS 1/4"
- Focus range: 8mm ~ 450mm
- Lens: high definition microscopy lens
- Dimension: 36mm (∅) x 142mm (length)
- PC interface: USB 2.0
- Mini USB cable: 1.2m
- Power consumption: max 2.5W
- Light source: 8 adjustable white LEDs

APPLICATIONS

- Dermatology
- Trichology
- Emergency Room
- Research & Scientific Analysis
- Telemedicine
- Telehealth
- Education



www.visio-tek.net For more information: info@visio-tek.net

TECHNICAL DATA SHEET